


VisionAire™

SERIES CONCENTRATORS



For a Look Into the
**Future of In-Home
Oxygen**



VISIONAIRE™ • VISIONAIRE™ 3
COMPACT OXYGEN CONCENTRATORS

AIRSEP®



A VisionAire-Series Concentrator for Non-Delivery Systems**

Choosing a model from one of the VisionAire-series compact concentrators—paired with AirSep's FreeStyle, the world's lightest and only wearable up to a 3-setting portable oxygen concentrator—packs a powerful scenario in a non-delivery system that is the most financially prudent in considering acquisition costs as well as economics for daily use. When the majority of hours of therapeutic use can be placed on this rugged, in-home stationary unit, the patient's POC is able to be reserved solely for out-of-home use, such as automobile and even airline travel.***

**See AirSep's Ultimate O₂ Combo brochure.

***FreeStyle has received the U.S. Federal Aviation Administration's (FAA) acceptance for onboard, in-flight use by oxygen passengers on commercial airlines via a 2006 amendment to SFAR 106. For a list of airlines that accept passengers with FreeStyle portable oxygen concentrators (POCs), reference www.airsep.com. Please check directly with the individual airlines for the most up-to-date information on their specific POC policies.

VisionAire

SPECIFICATIONS

Oxygen Concentration:	1-5 liters per minute at 90% + 5.5/-3% (Based on 70°F [21°C] at sea level)
Dimensions:	20.8 in. high, 14.1 in. wide, 11.5 in. deep (52.8 cm high, 35.8 cm wide, 29.2 cm deep)
Weight:	30 lb. shipping weight - 37 lb. (13.6 kg. shipping weight - 16.7 kg)
Electrical:	115 VAC, 60 Hz, 3 amps, 290 watts 220-240 VAC, 50 Hz, 1.5 amps, 290 watts 230 VAC, 60 Hz, 1.5 amps, 290 watts
Sound Level:	40 dBA
Alarms:	Power Failure High and Low Pressure High Temperature No Flow Low Oxygen Concentration (with optional Oxygen Monitor)
Relative Humidity:	Up to 95% (non-condensing)
WARRANTY:	Five-year warranty on parts. Other warranty options available.

INTERNATIONAL APPROVALS

 CE 0459
AirSep Corporation is ISO 9001:2008 and ISO 13485:2003-certified

ORDERING INFORMATION

AS098-1	VisionAire Oxygen Concentrator (115 V, 60 Hz unit)
AS098-2	VisionAire Oxygen Concentrator (220-240 V, 50 Hz unit)
AS098-3	VisionAire Oxygen Concentrator (230 V, 60 Hz unit)
AS098-4	VisionAire Oxygen Concentrator with oxygen monitor (115 V, 60 Hz unit)
AS098-5	VisionAire Oxygen Concentrator with oxygen monitor (220-240 V, 50 Hz unit)
AS098-6	VisionAire Oxygen Concentrator with oxygen monitor (230 V, 60 Hz unit)
AS098-707	VisionAire Oxygen Concentrator with air outlet (115 V, 60 Hz unit)
AS098-708	VisionAire Oxygen Concentrator with air outlet (220-240 V, 50 Hz unit)
AS098-709	VisionAire Oxygen Concentrator with air outlet (230 V, 60 Hz unit)
AS098-710	VisionAire Oxygen Concentrator with oxygen monitor and air outlet (115 V, 60 Hz unit)
AS098-711	VisionAire Oxygen Concentrator with oxygen monitor and air outlet (220-240 V, 50 Hz unit)
AS098-712	VisionAire Oxygen Concentrator with oxygen monitor and air outlet (230 V, 60 Hz unit)


Other model units available.

VisionAire 3

SPECIFICATIONS

Oxygen Concentration:	1-3 liters per minute at 90% + 5.5/-3% (Based on 70°F [21°C] at sea level)
Dimensions:	20.8 in. high, 14.1 in. wide, 11.5 in. deep (52.8 cm high, 35.8 cm wide, 29.2 cm deep)
Weight:	30 lb. shipping weight - 37 lb. (13.6 kg. shipping weight - 16.7 kg)
Electrical:	100 VAC, 50/60 Hz, 175 watts 115 VAC, 60 Hz, 175 watts 230 VAC, 50 Hz, 175 watts 230 VAC, 60 Hz, 175 watts
Sound Level:	<40 dBA
Alarms:	Power Failure High and Low Pressure High Temperature No Flow Low Oxygen Concentration (with optional Oxygen Monitor)
Relative Humidity:	Up to 95% (non-condensing)
WARRANTY:	Five-year warranty on parts. Other warranty options available.

INTERNATIONAL APPROVALS

 CE 0459
AirSep Corporation is ISO 9001:2008 and ISO 13485:2003-certified

ORDERING INFORMATION

AS096-1	VisionAire 3 Oxygen Concentrator (230 V, 50 Hz unit)
AS096-2	VisionAire 3 Oxygen Concentrator with oxygen monitor (230 V, 50 Hz unit)
AS096-3	VisionAire 3 Oxygen Concentrator with oxygen monitor, UK-style power cord (230 V, 60 Hz unit)
AS096-4	VisionAire 3 Oxygen Concentrator with oxygen monitor, Japan-style power cord (100 V, 50/60 Hz unit)
AS096-5	VisionAire 3 Oxygen Concentrator with oxygen monitor (230 V, 60 Hz unit)
AS096-6	VisionAire 3 Oxygen Concentrator (115 V, 60 Hz unit)
AS096-7	VisionAire 3 Oxygen Concentrator with oxygen monitor (115 V, 60 Hz unit)

Other model units available.



AirSep Corporation
Medical Products Division

401 Creekside Drive • Buffalo, New York 14228-2085 USA

Tel: 716.691.0202 • Fax: 716.691.4141 • E-mail: mpd@airsep.com • www.airsep.com